




Image Processing

Segmentation


- Segmentation is the **division** of an image into **discrete regions**.

Segmentation

Input




Semantic



Background

Instance



Background

Nucleus 1

Nucleus 2

Nucleus 3

...


How do we get segments?

- Thresholding-based
- Interactive tools based on classic machine learning
- Deep-learning based (Stardist, Cellpose)


Thresholding

The easiest way to segment an image is often by applying a *global threshold*.


This identifies pixels that are above or below a fixed threshold value, giving a *binary image as the output*.




Thresholding



Original, 8 bit grayscale
Blobs: Fiji example



Thresholded




Instance Segmented

Instance Segmentation in FIJI: keeping white (*connected-*) objects.

Thresholding method

Select only a **range** of **digital values** in the image.





Thresholding method

Select only a **range** of **digital values** in the image.

in **Fiji**: **Image > Adjust > Threshold...**

(cmd) + shift + t




The result of the thresholding process is a **Binary Mask**.



Thresholding method

Generate a **binary mask**.



Binary because the image has only **two** pixel values, **one** for the **selected pixels** and **one** for the “**discarded**” pixels.

In Fiji the two pixel values are **0** and **255**.



Segmentation with thresholding—exercises




7.1 DAPI segmentation with thresholding

skip “Analyze Particles”

Thresholding method


What can go wrong?

Usually, if you apply **thresholding** to the “**ORIGINAL**” image (the one you get out of the microscope), you won't be able to precisely **select all/only the pixels** you are interested in.



- **Fluorescence label** (e.g. DAPI)
- **Background** (uneven illumination, out-of-focus light, aberration, ...)
- **Noise** (detector read noise, Poisson noise, ...)

Interactive tools based on classic machine learning



Standalone: Ilastik



<https://www.ilastik.org/>

For histopathology



<https://qupath.github.io/>

As a Fiji plugin




<https://imagej.net/plugins/tws/>


WEKA: Waikato Environment for Knowledge Analysis: collection of free machine learning and data analysis software developed by the University of Waikato, NZ







Trainable WEKA Segmentation (in ImageJ/Fiji)



Trainable WEKA Segmentation (in ImageJ/Fiji)



Tips



Segmentation with pixel based classifier—exercises



7.2 DAPI segmentation with Weka

Breakpoint at 12: Create results